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LIFTmee LIP & EYE - Your Beauty Expert at home

LIFTmee LIP&EYE works with Micro Vibrations, these impulses get you valuable beautifying botanicals deep into the skin. Simultaneously, the blood and lymph circulation is stimulated, to give the skin a fresh and youthful appearance. Facial mimic muscles and small dry patches, especially in the sensitive lip and eye area can be visibly smoothed. The vibrating impulses, relaxes the skin surface and contributes to a general sense of well being.

Through its handy size LIFTmee LIP & EYE is not only an ideal travel companion but can also be used at the office, in front of the TV or comfortably before going to bed.

Application:

After cleaning the LIFTmee, apply sufficient serum or concentrate onto the skin and work in with small circular movements. The silver surface on the handle has to be touched at all times

Technical Specifications:

Voltage:	1.5 V DC
Energy consumption:	60 MW
Continuous running time:	6 hours
Size:	115 (L) x 24 (W) x 22 (H)
Battery:	1.5V AAA alkaline battery
Current:	0.1mA

Notes:

Please read the instruction manual before use. Never clean or place the appliance under running water. Use a damp cloth for cleaning.